

Clipper Vac® Attachment

Compatible with Wahl® Chromado Lithium



Step 1 Remove the blade from your clipper to have an open view of the upper housing.

Step 2 Clean both the top of the clipper and the attachment with alcohol to remove any residual oils. Open the small tube of silicone supplied with attachment and place a liberal amount on the entire bottom of the bottom of the black attachment.



Step 3 Place the Clipper Vac® Attachment, silicone side down, on the top of your clipper making sure to align the center mark on the attachment with the center mark on the body of the clipper. Press firmly to secure.

Step 4 Excess silicone should squeeze from under the attachment. Trim with a razor blade or other sharp object after it dries.

Note: *Allow 24 hours before using the clipper to ensure a proper seal on the clipper.*

M.D.C. Romani, Inc./Clipper Vac® is not affiliated with, authorized by or endorsed in any way by Wahl®.